

ICS Teaching Module: Electromyography in the assessment and therapy of lower urinary tract dysfunction in adults

J. Krhut, Ostrava, P.F.W.M. Rosier, B. Shelly, R.
Zachoval, P. Zvara



ICS

International Continence Society
Teaching Module

Electromyography in the assessment and therapy of lower urinary tract dysfunction in adults

Dr. Ömer Gülpınar

Principle of EMG

- Recording of the electrical activity from (striated) muscle with electrodes, to unveil function and innervation.

2 methods:

Needle EMG > Needle electrode(s): Inside muscle – motor unit

- Positive: - allows assessment of single action potentials
- Negative: - invasive - complex expertise in EMG required

Surface EMG > Surface (patch) electrode(s): On muscle – ‘whole’ muscle

- Positive: - non-invasive, less time- and money-consuming
- Negative: - less specific, less ‘detail’ - does not allow assessment of single action potentials

EMG tests in adult urology

- Needle EMG of external anal sphincter (EAS)
- Needle EMG of external urethral sphincter (EUS)
 - Monopolar
 - Bipolar (concentric)
 - Wire(s)
- Surface EMG of external anal sphincter (EAS)
- Surface EMG and sacral reflexes conductivity testing
- Surface EMG and biofeedback
- Surface EMG with cystometry and pressure/flow
- Surface EMG = 'kinesiographical EMG': with pair (or array) of electrodes over muscle

Needle EMG of external anal sphincter (EAS)

Principle:


- Recording of electrical activity of EAS
- Elements of muscle activity of the pelvic floor

Technique:

- Patient in lateral decubitus or lithotomy position
- EAS: Needle electrodes inserted bilaterally, approximately 0.5 cm lateral to the anus
- Assessment during maximal relaxation, during slight pelvic floor contraction, during maximal voluntary contraction and or during artificial bladder filling

Needle EMG of external anal sphincter (EAS)

- Evidence:
- Potentially useful to detect disturbances in neuroregulation of the pelvic floor muscles in patients with
 - lower motor neuron lesion
 - demyelinating diseases
 - with Parkinson disease
 - Multiple System Atrophy



Sakakibara R, et al.: J Neurol Neurosurg Psychiatry 2000; 68:25

Needle EMG of external urethral sphincter (EUS)

Principle:

- Direct recording of electrical activity of EUS

Technique:

- Patient in lateral decubitus or lithotomy position
- Needle electrodes inserted transperineally (♂) or transvaginally (♀) or transurethrally – via catheter
- Assessment during maximal relaxation, during slight pelvic floor contraction, during maximal voluntary contraction and or during artificial bladder filling

Needle EMG of external urethral sphincter (EUS)

Evidence:

- Limited evidence for role in clinical setting for EUS –EMG
- Some role in Fowler's (♀retention) syndrome
- Potentially useful in direct detection of electrical activity while bladder filling



Surface EMG of external anal sphincter (EAS)

Principle:

- Recording of muscle activity using surface (patch) electrodes or electrodes on cone or plug

Technique:

- Degreasing of the perianal skin
- 2 'active' electrodes adjusted bilaterally to the muco-cutaneous line + ground electrode
- Assessment of activity rest vs. contraction

Evidence:

- Tool to detect pelvic floor muscle activity or relaxation

Surface EMG and Sacral reflexes conductivity testing

Principle:


- Stimulation of the pudendal nerve to induce pelvic floor contraction to evaluate of bulbocavernosus (cliteroanal) reflex

Technique:

- Stimulation using electrode dorsal at the base on the penis (♂) or on clitoris (♀)
- The response recorded with surface or needle electrodes from the region of anal sphincter or bulbocavernosus muscle

Evidence:

- Absence or delay in response, suggest lower motor neuron impairment
- No relevant recent study which could support the role of this examination in the daily clinical work-up was found




Surface EMG and biofeedback

Principle:

- Detect the pelvic floor muscle activity and transform it into a visual and/or acoustic display in order convey the information to the patient

Technique:

- Surface electrodes are placed close to the anal sphincter or on an anal or intravaginal plug
- Recorded signal transformed into apparent sound or visual clue
- Allows the patient to better understand the functional status of the pelvic floor



Surface EMG and biofeedback

Clinical observations:

- Baseline between contractions – inconsistent and elevating
- Resting baseline – varies widely from session to session, especially when pain exists
- Erratic tracing without artifact or noise
- Patient has symptoms of overactive PFM – obstructed urination, defecation, pain
- Return to baseline after startle or frightening – overactive PFM is slow
- 2/3 of dysfunctional muscles will have normal resting baseline

Surface EMG and biofeedback

Evidence:

- Potentially useful for conservative treatment (PFM training) of stress urinary incontinence and OAB
- Little evidence regarding the use of EMG biofeedback as tool to help relax the pelvic floor muscles during micturition in adults

Surface EMG with cystometry

Principle:

- Recording of pelvic floor muscle activity during filling of the bladder.

Technique:

- Surface EMG (as on earlier slides)

Evidence:

- Introduced on the basis of expert opinion/ plausibility
- No (comparative) evidence
- Surface EMG may fail to reflect urethral (continence) function

Surface EMG with pressure flow studies

Principle and technique:


- Identical to surface EMG

Evidence:

- Introduced on the basis of expert opinion/ plausibility
- No (comparative) evidence
- Expert series demonstrating plausible results
- However: Large (n= 655) prospective cohort ♀ with EMG revised:
 - Many (51%) were not interpretable (but also)...
 - ...surface EMG failed to reflect EUS relaxation.

Surface EMG (adult)

- May be not interpretable - (technical) artefacts
- May give not plausible results - Not reflect relevant (EUS) activity - Smooth flow-rate and bizarre EMG



sEMG with urodynamic tests (adult)

Lacking practice standards for:

- Display: envelope; linear envelope; full wave; half wave
- Time scale
- Sampling: ... Hz; Filtering: moving average; root mean square
- Placing of active electrodes (♀; ♂)
- Impedance check(s) (cleaning of skin): <5 (or <10) kΩ
- Reference electrode - neutral! (Not on another muscle); trochanter; pelvic rim; sacrum
- Technical and clinical quality checks
- Analysis, interpretation and reporting

sEMG with urodynamic tests (adult)

- Not very invasive
- Not very time and cost consuming
- Without standard
- Without certainty of relevance
- May add confusion if artefacts are not acknowledged
- May be of help in (pelvic) muscle strength and control training and learning to relax

EMG tests in adult urology

- The concept of use of EMG methods in functional urology/urogynaecology and physiotherapy is supported by good theoretical basis
- Current value of EMG in diagnosis is however limited
- Currently EMG practice can only rarely play a decision making role in diagnostics of LUT partially due to lack of standards